Audéo Lumity Housing Colours

Slim Lumity **Housing Colours**









ActiveVent Receiver (MAV)
























The world's first intelligent hearing aid receiver



M Receiver Gain 51(2cc) • MPO 114(2cc)





Instant Fit - Domes























UP Receiver*





Acoustic Coupling



Mild to moderately-severe hearing loss

P Receiver **MAV** Receiver Gain 59(2cc) • MPO 122(2cc) Gain 50(2cc) • MPO 119(2cc)

Gain 71(2cc) • MPO 132(2cc) Moderately severe to profound hearing loss Moderate to severe hearing loss

More Phonak technology to enhance your experience





Marvel







CROS Paradise












10 NW 0





sonova



- Table top, Lanyard and Pointing modes
- NEW Headset mode for online calls
- NEW Volume indication on docking station





Table Top, Lanyard

& Select mode

Roger X Receiver Wirelessly install Roger X into Lumity hearing aids for RogerDirect





Features at a Glance

May 2023

Phonak Lumity.

PHONAK

A Sonova brand

Phonak Lumity Models

Audéo™ L



Audéo L-R L90/L70/L50/L30 Rechargeable



Audéo L-RT L90/L70/L50/L30 Rechargeable with telecoil



L90/L70/L50/L30 Rechargeable Waterproof* *up to 50cm

Audéo L-RL



Slim L

Slim L-R L90/L70 Rechargeable Left & right models with 7 degree curvature

Rechargeable Charging Options





- Contact charging
- Magnetic retention
- 3 hours to fully charge





- Compatible with Audéo L-R and Audéo L-RT
 - Induction charging
 - Built-in battery for 3 full charges on the go

Charger Case Go

Compatible with

Audéo L-RL



Audéo L-RL

3 hours to fully charge

Slim Charger Life Charger

- Compatible with Compatible with Slim L-R
- Compact, sleek design Induction Charging
 - Magnetic retention
 - 3 hours to fully charge

Wireless Accessories



Phonak Remote Control


the easy to use remote control



Phonak PartnerMic™

for better one to one conversations






into wireless headphones

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. 028-6697-02/2023-02/ Printed in Canada © Phonak Canada All rights reserved

SmartSpeech Technology

A collection of features working together to focus on improving speech understanding in all envrionments:

- Understand more speech in quiet and noise
- Spend less effort on listening
- Fully participate in challenging group conversations
- Enjoy hearing aids that automatically adjust



Features that shine in all technology levels

6


remium



adjustments via the

Advanced







Standard





FEATUREC DENIFFIE

FEATURES	BENEFITS
ActiveVent Performance	Intelligently balances hearing performance and comfort of own voice
Roger Technology	Provides significant improvement in speech understanding in noise and far distances
UltraZoom	In noise, zooms in on the voices you are facing
Motion Sensor Hearing	Automatically adjusts to optimize hearing while walking
Dynamic Noise Cancellation	Reduces background noise around you for improved hearing
StereoZoom 2.0	Seamlessly adjusts for more speech focus in the front in challenging noise
Speech Enhancer	Boosts soft and distant speech in quiet environments
SpeechSensor	Improves hearing from the side and behind when in noise
Calm situation	Provides listening comfort and clarity in quiet
Speech in Noise	Adapts to optimize speech in noisy environments
Comfort in Noise	Listening comfort in loud situations when no speech is present
Music	Offers a full and rich sound quality to your favorite tunes
Speech in Loud Noise	Uses StereoZoom 2.0 to help focus on speech in a loud noise environment
Speech in Car	Reduces car and road noise
Comfort in Echo	Reduces disturbing echo so you can hear speech clearly
Speech in 360° (manual)	Tracks and focuses on the dominant speaker from different sides
WindBlock	Suppresses wind noise so you can enjoy time outdoors
WhistleBlock	Eliminates whistling from the hearing aids before it is audible
NoiseBlock	Reduces background noise without reducing speech
SoundRelax	Cushions sudden loud or unpleasant impulse sounds
EchoBlock	Removes echo/reverberation to provide a cleaner sound
Tinnitus Balance	Noise generator to reduce the impact of your tinnitus
Tap Control	Pause/resume streaming, accept/end calls, and access voice assistant

FEATURES BENEFITS

	FEATURES	DENEFILS
SmartSpeech	ActiveVent Performance	Intelligently balances hearing performance and comfort of own voice
	Roger Technology	Provides significant improvement in speech understanding in noise and far distances
	UltraZoom	In noise, zooms in on the voices you are facing
	Motion Sensor Hearing	Automatically adjusts to optimize hearing while walking
	Dynamic Noise Cancellation	Reduces background noise around you for improved hearing
	StereoZoom 2.0	Seamlessly adjusts for more speech focus in the front in challenging noise
AutoSense OS 5.0	Calm situation	Provides listening comfort and clarity in quiet
	Speech in Noise	Adapts to optimize speech in noisy environments
	Comfort in Noise	Listening comfort in loud situations when no speech is present
	Music	Offers a full and rich sound quality to your favorite tunes
-	Speech in 360° (manual)	Tracks and focuses on the dominant speaker from different sides
	Speech in Loud Noise (manual)	Uses StereoZoom 2.0 to help focus on speech in a loud noise environment
	WindBlock	Suppresses wind noise so you can enjoy time outdoors
	WhistleBlock	Eliminates whistling from the hearing aids before it is audible
	NoiseBlock	Reduces background noise without reducing speech
	SoundRelax	Cushions sudden loud or unpleasant impulse sounds
	Tinnitus Balance	Noise generator to reduce the impact of your tinnitus
	Tap Control	Pause/resume streaming, accept/end calls, and access voice assistant

FEATURES BENEFITS

ActiveVent Performance	Intelligently balances hearing performance and comfort of own voice
Roger Technology	Provides significant improvement in speech understanding in noise and far distances
UltraZoom	In noise, zooms in on the voices you are facing
Motion Sensor Hearing	Automatically adjusts to optimize hearing while walking
Calm situation	Provides listening comfort and clarity in quiet
Speech in Noise	Adapts to optimize speech in noisy environments
Comfort in Noise	Listening comfort in loud situations when no speech is present
Music	Offers a full and rich sound quality to your favorite tunes
WindBlock	Suppresses wind noise so you can enjoy time outdoors
WhistleBlock	Eliminates whistling from the hearing aids before it is audible
NoiseBlock	Reduces background noise without reducing speech
Tinnitus Balance	Noise generator to reduce the impact of your tinnitus

Latzel, M., Lesimple, C., & Woodward, J. (2023). Speech Enhancer significantly reduces listening effort and increases intelligibility for speech at a distance. Phonak Field Study News. Retrieved from www.phonak.com/evidence.